

The Fully Adjustable Precision Torque Wrench

The Precise™ Adjustable Torque Wrench



Features

- CE Mark
- Autoclavable
- Ratchet Mechanism
-

9 Torque Settings

15, 20, 25, 30, 35, 40, 45, 50 and 60 Ncm

- Torque settings lock into select Location for accuracy and repeatability
- Ratchet clutch disengages (slips) when torque is applied to avoid over torquing
- Torque settings are provided for use with all implant designs
- * Wrench head accepts the standard 4mm square connection



IMPLANT SYSTEM	DRIVERS																	
	.035 HCT	.035 HCLT	.040 HSM	.040 HSM L	.048 HNT	.048 HNL T	.050 HTT	.050 HTL T	.062 HS	.062 HSL	.080 HJ	.080 HJ L	.080 HFD	.177 HTSI	.050 SQ CER	TORX PROC	TORX ITI	TGD
NobelBiocare	•	•			•	•							•		•	•		
Keystone Dental	•	•			•	•			•	•					•			
Steri-Oss							•	•										
Straumann/ITI																	•	•
Friadent	•	•	•	•	•	•												
Sterngold / Impla-Med					•	•							•					
Implant Innovations/31	•	•			•	•							•		•			•
Implant Direct							•	•				•	•					
Zimmer									•	•								
BioHorizons	•	•					•	•										
Astra							•	•										
Camlog							•	•										